Paths of analysis*

Synthia

March 3, 2022

1 Analysis parameters

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: FGI, FGI with protections

Max. paths returned: 5

Max. iterations: 300

Commercial:

1. Max. molecular weight - 1000 g/mol

2. Max. price - 1000 \$/g

Published:

- 1. Max. molecular weight 1000 g/mol
- 2. Popularity 10

My Stockroom:

1. Max. molecular weight - 1000 g/mol

Reaction scoring formula: TUNNEL_COEF*FGI_COEF*STEP*20+1000 000*(CONFLICT+NON SELECTIVITY+FILTERS+PROTECT)

Chemical scoring formula: SMALLER^ 3,SMALLER^ 1.5

Min. search width: 400

Max. reactions per product: 60

Strategies: none selected

^{*}The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

FGI Coeff: 0

JSON Parameters: {}

2 Paths

3 paths found. Paths are sorted by score. Reactions are sorted in appearance order for each path.

2.1 Path 1

Score: 576.15




Figure 1: Outline of path 1

2.1.1 Photoredox Cross-Electrophile Coupling of Unactivated Alkyl Bromides

$$+ H_2N$$

$$Br$$

$$S + S$$

$$?$$

Substrates:

1. 2-Bromopropane - available at Sigma-Aldrich

2. 5-Bromo-2-isopropoxyaniline - available at Sigma-Aldrich

Products:

1. CC(C)Oc1ccc(C(C)C)cc1N

Typical conditions: [Ir]-photocat.[Ni]-cat.TTMSS.base.blue light

Protections: none


Yield: good

Reference: 10.1021/jacs.6b04818 and 10.1016/j.bbrc.2020.04.028 and 10.1021/ac-

smedchem lett. 8b00183

Retrosynthesis ID: 31016940

2.1.2 Bromination of aromatic compounds



Substrates:

1. CC(C)Oc1ccc(C(C)C)cc1N

Products:

1. CC(C)Oc1cc(Br)c(C(C)C)cc1N

Typical conditions: Br2.Fe

Protections: none

Yield: good

Reference: 10.1021/acs.accounts.6b00120

2.1.3 Photoredox Cross-Electrophile Coupling of Unactivated Alkyl Bromides

Substrates:

- 1. CC(C)Oc1cc(Br)c(C(C)C)cc1N

Products:

1. CC(C)Oc1cc(C2CC2)c(C(C)C)cc1N

Typical conditions: [Ir]-photocat.[Ni]-cat.TTMSS.base.blue light

Protections: none

Yield: good

Reference: 10.1021/jacs.6b04818 and 10.1016/j.bbrc.2020.04.028 and 10.1021/ac-j.bbrc.2020.04.028 and 10.1021/ac-j.bbrc.2020.04.028

smedchemlett.8b00183

Retrosynthesis ID: 31016940

2.1.4 Amination of aryl bromides

Substrates:

1. 2-Bromo-N,N-dimethylbenzenesulphonamide 1g pack - Combi-Blocks

2. Benzylamine - available at Sigma-Aldrich

Products:

1. CN(C)S(=O)(=O)c1ccccc1NCc1ccccc1

Typical conditions: Pd.ligand.base or CuI.ligand.base

Protections: none


Yield: good

Reference: 10.1021/ja903049z and 10.1021/jo060945k and 10.1021/jo060190h

and 10.1039/B923255A and 10.1021/jm8003625 and 10.1021/jo9006738

Retrosynthesis ID: 28544

2.1.5 Synthesis of haloarenes via triflates



Substrates:

1. 5-chloro-uracil - Combi-Blocks

Products:

 $1. \ O{=}c1[nH]c(F)ncc1Cl$

Typical conditions: 1.Tf2O 2. [Pd].MX

Protections: none
Yield: moderate

Reference: 10.1016/j.tetasy.2012.04.008 and WO2007/136577 (p46) and

10.1021/ol202098h and 10.1021/ol402859k and 10.1021/jacs.5b09308

2.1.6 Amination of pyridones

Substrates:

- 1. CN(C)S(=O)(=O)c1ccccc1NCc1ccccc1
- 2. O=c1[nH]c(F)ncc1Cl

Products:

 $1. \ \mathrm{CN(C)S(=O)(=O)c1ccccc1N(Cc1ccccc1)c1nc(F)ncc1Cl}$

Typical conditions: 1.PCl5.2.amine

Protections: none
Yield: moderate

Reference: 10.1021/jm300780p AND 10.3390/molecules170910902 AND

10.1021/jm00392a017

Retrosynthesis ID: 14895

2.1.7 Nucleophilic aromatic substitution

Substrates:

1. CN(C)S(=O)(=O)c1ccccc1N(Cc1ccccc1)c1nc(F)ncc1Cl

 $2. \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1N$

Products:

 $1. \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1\\$

Typical conditions: Solvent


Protections: none

Yield: good

Reference: 10.1002/9781118093559.ch4

Retrosynthesis ID: 49476

2.1.8 Debenzylation



Substrates:

 $1. \ \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1$

Products:

 $1. \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(Nc2cccc2S(=O)(=O)N(C)C)n1 \\$

Typical conditions: H2. Pd/C or Pd(OH)2

Protections: none

Yield: good

Reference: DOI: 10.1002/1521-3773(20020603)41:11<1895::AID-ANIE1895>3.0.CO;2-3 and 10.1021/jo400589j and 10.1021/jm8012932 (SI,page S6) and 10.1080/00397911.2016.1261164

Retrosynthesis ID: 9995661

2.2 Path 2

Score: 639.39





Figure 2: Outline of path 2

2.2.1 Alkylation of amines with alkyl chlorides (PTC conditions)



Substrates:

- 1. a-Chlorotoluene available at Sigma-Aldrich
- 2. 2-Amino-N,N-dimethylbenzenesulfonamide Combi-Blocks

Products:

1. CN(C)S(=O)(=O)c1ccccc1NCc1ccccc1

Typical conditions: NaOH.water.PTC-catalyst


Protections: none
Yield: moderate

 $\textbf{Reference:} \quad 10.1080/00397911.2013.828077 \quad \text{and} \quad 10.1002/ejoc.201200202 \quad \text{and} \quad 10.100200202 \quad \text{and} \quad 10.100200202020202 \quad \text{and} \quad 10.100200202020202020202020202$

10.1080/10799893.2019.1585453 and 10.1248/cpb.c14-00754

Retrosynthesis ID: 4785

2.2.2 Synthesis of halopyrimidines



Substrates:

1. 5-Chlorouracil - available at Sigma-Aldrich

Products:

1. O=c1[nH]c(Br)ncc1Cl

Typical conditions: POCl3.or.POBr3 or Ph3P.NBS/NCS

Protections: none
Yield: moderate

Reference: 10.1002/1522-2675(20010516)84:5<1112::AID-HLCA1112>3.0.CO;2-8 and 10.1039/c5dt02185e and 10.1016/j.molstruc.2012.07.010 and 10.1021/jo01293a022 and 10.1002/ejoc.200400209

2.2.3 Photoredox Cross-Electrophile Coupling of Unactivated Alkyl Bromides

$$Br \longrightarrow H_2N \longrightarrow H_2N \longrightarrow Br$$

$$S + S \longrightarrow P$$

Substrates:

1. 2-Bromopropane - available at Sigma-Aldrich

2. 5-Bromo-2-isopropoxyaniline - available at Sigma-Aldrich

Products:

1. CC(C)Oc1ccc(C(C)C)cc1N

Typical conditions: [Ir]-photocat.[Ni]-cat.TTMSS.base.blue light

 ${\bf Protections:}\ {\bf none}$

Yield: good

Reference: 10.1021/jacs.6b04818 and 10.1016/j.bbrc.2020.04.028 and 10.1021/ac-

smedchem lett. 8b00183

Retrosynthesis ID: 31016940

2.2.4 Amination of pyridones

 ${\bf Substrates:}$

1. CN(C)S(=O)(=O)c1ccccc1NCc1ccccc1

2. O=c1[nH]c(Br)ncc1Cl

Products:

1. CN(C)S(=O)(=O)c1ccccc1N(Cc1ccccc1)c1nc(Br)ncc1Cl

Typical conditions: 1.PCl5.2.amine


Protections: none
Yield: moderate

Reference: 10.1021/jm300780p AND 10.3390/molecules170910902 AND

10.1021/jm00392a017

Retrosynthesis ID: 14895

2.2.5 Iodination of aromatic compounds



Substrates:

1. CC(C)Oc1ccc(C(C)C)cc1N

Products:

1. CC(C)Oc1cc(I)c(C(C)C)cc1N

Typical conditions: I2 or other iodinating agent e.g. NIS

Protections: none

Yield: good

Reference: DOI: 10.1039/C5SC00964B and 10.1016/j.tetlet.2005.05.117 and 10.1007/s11178-005-0256-1

Retrosynthesis ID: 10697

2.2.6 One pot synthesis of aryl-alkyl sulfides

Substrates:

1. Hydrogen sulfide - available at Sigma-Aldrich

2. Methyl bromide - available at Sigma-Aldrich

3. CN(C)S(=O)(=O)c1ccccc1N(Cc1ccccc1)c1nc(Br)ncc1Cl

Products:

1. CSc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1

Typical conditions: nBuLi.THF.-78C.then.S.then AlkBr.to.rt

Protections: none

Yield: good

Reference: 10.1021/jo049758h Retrosynthesis ID: 5320

2.2.7 Suzuki coupling of cyclopropanotrifluoroboranes and aryl iodides



Substrates:

1. CC(C)Oc1cc(I)c(C(C)C)cc1N

2. Potassium cyclopropyltrifluoroborate - available at Sigma-Aldrich

Products:

1. CC(C)Oc1cc(C2CC2)c(C(C)C)cc1N

Typical conditions: K3PO4.Pd(OAc)2.ruphos.H2O.toluene.110C

Protections: none

Yield: good

Reference: 10.1002/ejoc.201100119 and WO2014006066 p.26 and

WO2014184275 p.107

Retrosynthesis ID: 31016521

2.2.8 Substitution of 2-thiomethylpyrimidines with amines

Substrates:

1. CC(C)Oc1cc(C2CC2)c(C(C)C)cc1N

 $2. \ CSc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1 \\$

Products:

 $1. \ \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1$

Typical conditions: K2CO3.DMF

Protections: none

Yield: good

Reference: 10.1021/jm980222w AND 10.1016/j.cclet.2014.10.007 AND

10.1002/jhet.5570280520 AND 10.1080/00397910701396930

Retrosynthesis ID: 14935

2.2.9 Debenzylation

Substrates:

 $1. \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1$

Products:

 $1. \ \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(Nc2cccc2S(=O)(=O)N(C)C)n1$

Typical conditions: H2. Pd/C or Pd(OH)2

Protections: none

Yield: good

Reference: DOI: 10.1002/1521-3773(20020603)41:11<1895::AID-ANIE1895>3.0.CO;2-3 and 10.1021/jo400589j and 10.1021/jm8012932 (SI,page S6) and 10.1080/00397911.2016.1261164

Retrosynthesis ID: 9995661

2.3 Path 3

Score: 741.80




Figure 3: Outline of path 3

2.3.1 Photoredox Cross-Electrophile Coupling of Unactivated Alkyl Bromides

$$Br \longrightarrow H_2N \longrightarrow H_2N \longrightarrow P$$

Substrates:

- 1. 2-Bromopropane available at Sigma-Aldrich
- 2. 5-Bromo-2-isopropoxyaniline available at Sigma-Aldrich

Products:

1. CC(C)Oc1ccc(C(C)C)cc1N

 $\textbf{Typical conditions:} \ [Ir]-photocat.[Ni]-cat.TTMSS.base.blue \ light$

Protections: none

Yield: good

Reference: 10.1021/jacs.6b04818 and 10.1016/j.bbrc.2020.04.028 and 10.1021/ac-j.bbrc.2020.04.028 and 10.1021/ac-j.bbrc.2020.04.028

smedchem lett. 8b00183

Retrosynthesis ID: 31016940

2.3.2 Chlorination of aromatic compounds

$$H_2N$$
 H_2N
 CI

${\bf Substrates:}$

1. CC(C)Oc1ccc(C(C)C)cc1N

Products:

1. CC(C)Oc1cc(Cl)c(C(C)C)cc1N

Typical conditions: Cl2 or other chlorinating agent like NCS

Protections: none
Yield: moderate

Reference: DOI: 10.1007/s11178-005-0256-1

${\bf 2.3.3} \quad {\bf Suzuki\ coupling\ of\ cyclopropanotrifluor oboranes\ and\ aryl\ chlorides}$

Substrates:

- 1. CC(C)Oc1cc(Cl)c(C(C)C)cc1N
- 2. Potassium cyclopropyltrifluoroborate available at Sigma-Aldrich

Products:

1. CC(C)Oc1cc(C2CC2)c(C(C)C)cc1N

Typical conditions: K2CO3.Pd(OAc)2.H2O.100C

Protections: none

Yield: good

Reference: 10.1021/jo801269m and WO2012101066 p.177 and WO2011/55115

p.72

Retrosynthesis ID: 31016519

2.3.4 Alkylation of amines with alkyl chlorides (PTC conditions)

CI
$$+ CI$$

$$NH_2$$

$$N \rightarrow CI$$

$$S + S \rightarrow 2$$

Substrates:

1. 2,5-Dichloropyrimidin-4-amine - Combi-Blocks

2. a-Chlorotoluene - available at Sigma-Aldrich

Products:

1. Clc1ncc(Cl)c(NCc2cccc2)n1

Typical conditions: NaOH.water.PTC-catalyst


Protections: none
Yield: moderate

Reference: 10.1080/00397911.2013.828077 and 10.1002/ejoc.201200202 and

10.1080/10799893.2019.1585453 and 10.1248/cpb.c14-00754

Retrosynthesis ID: 4785

2.3.5 Nucleophilic aromatic substitution



Substrates:

1. Clc1ncc(Cl)c(NCc2cccc2)n1

Products:

1. Fc1ncc(Cl)c(NCc2cccc2)n1

Typical conditions: CsF.DMF

Protections: none
Yield: moderate

Reference: 10.1016/j.tetlet.2015.09.057 and 10.1016/j.cattod.2012.02.063 and

10.1016/j.ejmech.2015.01.034

2.3.6 Chan-Lam Coupling

Substrates:

- $1. \ \, Fc1ncc(Cl)c(NCc2cccc2)n1$
- 2. 2-(N,N-DIMETHYLSULPHAMOYL)BENZENEBORONIC ACID Combi-Blocks

Products:

1. CN(C)S(=O)(=O)c1ccccc1N(Cc1ccccc1)c1nc(F)ncc1Cl

 $\textbf{Typical conditions:} \ \mathrm{Cu(II)Ac.DCM}$

Protections: none

Yield: good

Reference: 10.5012/bkcs.2012.33.5.1785 and 10.1039/C8GC02611D and

10.1055/s-2008-1032184

Retrosynthesis ID: 31015959

2.3.7 Nucleophilic aromatic substitution

Substrates:

1. Methanethiol - available at Sigma-Aldrich

 $2. \ \ CN(C)S(=O)(=O)c1ccccc1N(Cc1ccccc1)c1nc(F)ncc1Cl$

Products:

1. CSc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1

Typical conditions: NaH.THF.0-80 C or K2CO3.DMF.110 C

Protections: none


Yield: good

10.1016/j.tetlet.2015.10.008 p. Reference: 6479,6483 and

10.1016/j.ejmech.2016.06.056 p. 82, 85

Retrosynthesis ID: 49475

2.3.8 Substituion of 2-thiomethylpyrimidines with amines



Substrates:

1. CC(C)Oc1cc(C2CC2)c(C(C)C)cc1N

2. CSc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1

Products:

 $1. \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1\\$

Typical conditions: K2CO3.DMF

Protections: none

Yield: good

Reference: 10.1021/jm980222w AND 10.1016/j.cclet.2014.10.007 AND

10.1002/jhet.5570280520 AND 10.1080/00397910701396930

Retrosynthesis ID: 14935

2.3.9 Debenzylation

Substrates:

 $1. \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(N(Cc2cccc2)c2cccc2S(=O)(=O)N(C)C)n1\\$

Products:

 $1. \ \ CC(C)Oc1cc(C2CC2)c(C(C)C)cc1Nc1ncc(Cl)c(Nc2cccc2S(=O)(=O)N(C)C)n1$

Typical conditions: H2. Pd/C or Pd(OH)2

Protections: none

Yield: good

Reference: DOI: 10.1002/1521-3773(20020603)41:11<1895::AID-ANIE1895>3.0.CO;2-3 and 10.1021/jo400589j and 10.1021/jm8012932 (SI,page SG) and 10.1020/00207011 2015 1001151

S6) and 10.1080/00397911.2016.1261164