Paths of analysis* Analysis 2

Synthia

March 3, 2022

Analysis parameters 1

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: FGI, FGI with protections

Max. paths returned: 5

Max. iterations: 300

Commercial:

- 1. Max. molecular weight 1000 g/mol
- 2. Max. price 1000 \$/g

Published:

- 1. Max. molecular weight 1000 g/mol
- 2. Popularity 10

My Stockroom:

1. Max. molecular weight - 1000 g/mol

Reaction scoring formula: TUNNEL COEF*FGI COEF*STEP*20+1000 000*(CONFLICT+NON SELECTIVITY+FILTERS+PROTECT)

Chemical scoring formula: SMALLER^ 3,SMALLER^ 1.5

Min. search width: 400

Max. reactions per product: 60

Strategies: none selected

^{*}The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

FGI Coeff: 0

JSON Parameters: {}

2 Paths

3 paths found. Paths are sorted by score. Reactions are sorted in appearance order for each path.

2.1 Path 1

Score: 77.96




Figure 1: Outline of path 1

2.1.1 Synthesis of N-arylamides from arenediazonium salts

$$H_2N \longrightarrow 0$$
 $H_2N \longrightarrow 0$
 $H_2N \longrightarrow 0$

Substrates:

- 1. 4-Bromoaniline available at Sigma-Aldrich
- 2. Calcium nitrite solution available at Sigma-Aldrich
- 3. 2-Chloronicotinamide Combi-Blocks

Products:

1. n-(4-bromo-phenyl)-2-chloro-nicotinamide - Enamine

Typical conditions: 1) HCl.NaNO2 2) CuI.TBAI.N,N'-dimethylethane-1,2-

 ${\it diamine.} K2CO3.DMSO.110C$

Protections: none
Yield: moderate

Reference: DOI: 10.1055/s-0034-1378556

Retrosynthesis ID: 1922

2.1.2 Pd-catalyzed coupling of aryl bromides and alkenes

Substrates:

1. n-(4-bromo-phenyl)-2-chloro-nicotinamide - Enamine

2. b-Methallyl chloride - available at Sigma-Aldrich

Products:

1. CC(C)(CCl)c1ccc(NC(=O)c2cccnc2Cl)cc1

Typical conditions: [Pd(P(o-tol)3OAc]2.NaOCHO.TBAB.TEA

Protections: none
Yield: moderate

Reference: DOI: 10.1038/nature05569

Retrosynthesis ID: 8840

2.1.3 Buchwald-Hartwig amination

Substrates:

1. 4-Picolylamine - available at Sigma-Aldrich

2. CC(C)(CCl)c1ccc(NC(=O)c2cccnc2Cl)cc1

Products:

 $1. \ CC(C)(CCl)c1ccc(NC(=O)c2cccnc2NCc2ccncc2)cc1 \\$

 ${\bf Typical\ conditions:}\ {\bf PdCl2.NaOtBu.dioxane.heat}$

Protections: none

Yield: good

Reference: 10.1021/acs.oprd.9b00161 and 10.1002/anie.201904795 and

10.1021/acs. chemrev. 6b00512

Retrosynthesis ID: 10319

2.2 Path 2

Score: 112.45




Figure 2: Outline of path 2

2.2.1 Amide coupling

Substrates:

- 1. 4-Bromoaniline available at Sigma-Aldrich
- 2. 2-Nitronicotinic acid Combi-Blocks

Products:

1. O=C(Nc1ccc(Br)cc1)c1cccnc1[N+](=O)[O-]

Typical conditions: DCC.DCM or EDC.DCM or SOCl2.DCM

Protections: none

 $\mathbf{Yield}: \mathbf{good}$

Reference: 10.1021/cr100048w and 10.1039/B701677H and 10.1039/C5RA24527C and 10.3727/0000000006783981206 and 10.1021/np060007f and 10.1021/j000012a058 and 10.1016/j.bmcl.2007.08.037 and 10.1039/C0OB00355G and 10.1021/jm500031w (p.3056) and 10.1016/j.tet.2011.03.046

2.2.2 Pd-catalyzed coupling of aryl bromides and alkenes

Substrates:

- 1. O=C(Nc1ccc(Br)cc1)c1cccnc1[N+](=O)[O-]
- 2. b-Methallyl chloride available at Sigma-Aldrich

Products:

 $1. \ \mathrm{CC(C)(CCl)c1ccc(NC(=O)c2cccnc2[N+](=O)[O-])cc1}$

 $\textbf{Typical conditions:} \ [Pd(P(o\text{-}tol)3OAc]2.NaOCHO.TBAB.TEA$

Protections: none
Yield: moderate

Reference: DOI: 10.1038/nature05569

2.2.3 Palladium-catalyzed reduction of nitro group

Substrates:

 $1. \ \mathrm{CC(C)(CCl)c1ccc(NC(=O)c2cccnc2[N+](=O)[O-])cc1}$

Products:

 $1. \ \mathrm{CC(C)(CCl)c1ccc(NC(=O)c2cccnc2N)cc1}$

Typical conditions: $\mathrm{H2.Pd/C}$

Protections: none
Yield: moderate

Reference: DOI: 10.1002/anie.200501738 and 10.1002/anie.200352175 and 10.1016/j.tetlet.2015.05.004 and 10.3390/molecules13061427 and 10.1016/S0968-0006(00)00450.00

0896(03)00459-0

2.2.4 Synthesis of secondary amines

Substrates:

- $1. \ \mathrm{CC(C)(CCl)c1ccc(NC(=O)c2cccnc2N)cc1}$
- 2. Isonicotinaldehyde available at Sigma-Aldrich

Products:

1. CC(C)(CCl)c1ccc(NC(=O)c2cccnc2NCc2ccncc2)cc1

 ${\bf Typical\ conditions:}\ {\bf sodium\ triacetoxyborohydride.} {\bf dichloromethane}$

Protections: none

Yield: good

Reference: DOI: 10.1021/ed077p270

Retrosynthesis ID: 245716

2.3 Path 3

Score: 119.74





Figure 3: Outline of path 3

2.3.1 Bromination of aromatic compounds



Substrates:

1. Neophyl chloride - available at Sigma-Aldrich

Products:

1. p-brom-(chlor-tert-butyl)-benzol


Typical conditions: Br2.Fe

Protections: none

 $\bf Yield: \ good$

Reference: 10.1021/acs.accounts.6b00120

2.3.2 Amination of aryl bromides



Substrates:

1. p-brom-(chlor-tert-butyl)-benzol

Products:

 $1. \ \, p\text{-}(1\text{-}chloro\text{-}2\text{-}methyl\text{-}2\text{-}propyl) aniline$

 ${\bf Typical\ conditions:}\ {\bf Pd.ligand.base\ or\ CuI.ligand.base}$

Protections: none

 $\mathbf{Yield}: \mathbf{good}$

Reference: 10.1021/ja903049z and 10.1021/jo060945k and 10.1021/jo060190h and 10.1039/B923255A and 10.1021/jm8003625 and 10.1021/jo9006738

2.3.3 Arylation of carbamates with aryl bromides

Substrates:

- 1. N-Boc-4-aminomethylpyridine available at Sigma-Aldrich
- 2. 2-Bromonicotinonitrile Combi-Blocks

Products:

1. CC(C)(C)OC(=O)N(Cc1ccncc1)c1ncccc1C#N

Typical conditions: Base.[Pd].catalyst.dioxane.heat or CuI.diamine.base.DMF.heat

Protections: none

 $\mathbf{Yield}: \mathbf{good}$

Reference: 10.1016/j.tetlet.2014.03.016 and 10.1021/ja012610k and

10.1021/ol016208m and 10.1021/ol502322c

2.3.4 Base hydrolysis of nitriles to carboxylic acids

Substrates:

1. CC(C)(C)OC(=O)N(Cc1ccncc1)c1ncccc1C#N

Products:

1. CC(C)(C)OC(=O)N(Cc1ccncc1)c1ncccc1C(=O)O

Typical conditions: NaOH.heating.H2O

Protections: none
Yield: moderate

Reference: 10.1002/1099-0690(200111)2001:22<4207::AID-EJOC4207>3.0.CO;2-

3 and 10.1021/acs.jmedchem.5b00702 and 10.1016/j.bmc.2011.07.045

2.3.5 Amide coupling

Substrates:

- 1. CC(C)(C)OC(=O)N(Cc1ccncc1)c1ncccc1C(=O)O
- 2. p-(1-chloro-2-methyl-2-propyl)aniline

Products:

 $1. \ CC(C)(C)OC(=O)N(Cc1ccncc1)c1ncccc1C(=O)Nc1ccc(C(C)(C)CC1)cc1$

Typical conditions: DCC.DCM or EDC.DCM or SOC12.DCM

Protections: none

Yield: good

Reference: 10.1021/cr100048w and 10.1039/B701677H and 10.1039/C5RA24527C and 10.3727/0000000006783981206 and 10.1021/np060007f and 10.1021/jo00012a058 and 10.1016/j.bmcl.2007.08.037 and 10.1039/C0OB00355G and 10.1021/jm500031w (p.3056) and 10.1016/j.tet.2011.03.046

2.3.6 Boc removal

Substrates:

 $1. \ CC(C)(C)OC(=O)N(Cc1ccncc1)c1ncccc1C(=O)Nc1ccc(C(C)(C)CC1)cc1$

Products:

 $1. \ \mathrm{CC(C)(CCl)c1ccc(NC(=O)c2cccnc2NCc2ccncc2)cc1}$

Typical conditions: TFA.DCM

Protections: none

Yield: good

Reference: 10.1016/j.bmc.2015.11.006 and 10.1021/jo047752m and

10.1016/j.tetlet.2007.09.003