# Paths of analysis\*

# Synthia

March 3, 2022

# 1 Analysis parameters

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: FGI, FGI with protections

Max. paths returned: 5

Max. iterations: 300

### Commercial:

1. Max. molecular weight - 1000 g/mol

2. Max. price - 1000 \$/g

#### Published:

- 1. Max. molecular weight 1000 g/mol
- 2. Popularity 10

### My Stockroom:

1. Max. molecular weight - 1000 g/mol

**Reaction scoring formula:** TUNNEL\_COEF\*FGI\_COEF\*STEP\*20+1000 000\*(CONFLICT+NON SELECTIVITY+FILTERS+PROTECT)

Chemical scoring formula: SMALLER^ 3,SMALLER^ 1.5

Min. search width: 400

Max. reactions per product: 60

Strategies: none selected

<sup>\*</sup>The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

### FGI Coeff: 0

JSON Parameters: {}

# 2 Paths

 $2~{\rm paths}$  found. Paths are sorted by score. Reactions are sorted in appearance order for each path.

# 2.1 Path 1

Score: 133.42





Figure 1: Outline of path 1

# 2.1.1 Balz-Schiemann Reaction



# Substrates:

 $1. \ \, \text{4-Bromopyrimidin-2-amine} \, \text{-} \quad \, \textit{Combi-Blocks}$ 

# **Products:**

1. 4-Bromo-2-fluoropyrimidine - Enamine

Typical conditions: NaNO2.HF-pyridine.-25 to 0C  $\,$ 

 ${\bf Protections:}\ {\bf none}$ 

Yield: moderate

**Reference:** 10.1021/jm100432w and 10.1021/jm5008177 and 10.1021/o1401540k and 10.1021/jm401551n (main text and SI, page S8) and 10.1515/chempap-2016-0033 and 10.1021/jm0311442 and 10.1016/j.jfluchem.2007.03.012 and 10.1021/jo00185a023

Retrosynthesis ID: 29906

# 2.1.2 Photoredox Cross-Electrophile Coupling of Unactivated Alkyl Bromides

### Substrates:

1. tert-Butyl bromide - available at Sigma-Aldrich

2. 3-Bromo-4-fluoroaniline - Combi-Blocks

#### **Products:**

1. CC(C)(C)c1cc(N)ccc1F

Typical conditions: [Ir]-photocat.[Ni]-cat.TTMSS.base.blue light

Protections: none

Yield: good

**Reference:** 10.1021/jacs.6b04818 and 10.1016/j.bbrc.2020.04.028 and 10.1021/ac-

smedchemlett.8b00183

# 2.1.3 Meerwein coupling of diazonium salt with heteroaryl

### Substrates:

- 1. 4-Bromo-2-fluoropyrimidine Enamine
- 2. CC(C)(C)c1cc(N)ccc1F

# Products:

1. CC(C)(C)c1cc(-c2nc(F)ncc2Br)ccc1F

Typical conditions: 1) HCl.NaNO2 2) [Ru(bpy)3Cl2]\*6H2O.45W

bulb.H2O.rt

 ${\bf Protections:}\ {\bf none}$ 

Yield: moderate

**Reference:** 10.1002/chem.201304120

Retrosynthesis ID: 10001815

# 2.1.4 Nucleophilic aromatic substitution

Substrates:

- 1. 5-((4-Ethylpiperazin-1-yl)methyl)pyridin-2-amine available at Sigma-
- 2. CC(C)(C)c1cc(-c2nc(F)ncc2Br)ccc1F

### Products:

 $1. \ \ CCN1CCN(Cc2ccc(Nc3ncc(Br)c(-c4ccc(F)c(C(C)(C)C)c4)n3)nc2)CC1$ 

Typical conditions: Solvent

Protections: none

Yield: good

Reference: 10.1002/9781118093559.ch4

Retrosynthesis ID: 49476

# 2.2 Path 2

Score: 185.81




Figure 2: Outline of path 2

# 2.2.1 Photoredox Cross-Electrophile Coupling of Unactivated Alkyl Bromides



#### Substrates:

1. tert-Butyl bromide - available at Sigma-Aldrich

2. 3-Bromo-4-fluoroaniline - Combi-Blocks

### **Products:**

1. CC(C)(C)c1cc(N)ccc1F

 $\textbf{Typical conditions:} \ [Ir]-photocat. [Ni]-cat. TTMSS. base. blue \ light$ 

Protections: none


Yield: good

**Reference:** 10.1021/jacs.6b04818 and 10.1016/j.bbrc.2020.04.028 and 10.1021/ac-j.bbrc.2020.04.028 and 10.1021/ac-j.bbrc.2020.04.028

smedchemlett.8b00183

Retrosynthesis ID: 31016943

# 2.2.2 Synthesis of boronic ester derivatives from amines



### Substrates:

1. CC(C)(C)c1cc(N)ccc1F

2. Bis(pinacolato)diboron - available at Sigma-Aldrich

# **Products:**

1. CC(C)(C)c1cc(B2OC(C)(C)C(C)(C)O2)ccc1F

Typical conditions: eg. HCl. t-BuONO. 0C or photoredox condition

Protections: none

Yield: good

**Reference:** 10.1021/acs.accounts.7b00566 and 10.1016/j.tetlet.2016.02.097 and 10.1002/anie.200905824 and 10.3184/174751918X15362432558247 and 10.1055/s-0031-1290960 and 10.1016/j.tetlet.2017.08.060 and 10.1021/jo3018878

Retrosynthesis ID: 2205

# 2.2.3 Suzuki coupling of arylboronic acids pinacol esters with aryl bromides

#### Substrates:

1. 4-Bromopyrimidin-2-amine - Combi-Blocks

2. CC(C)(C)c1cc(B2OC(C)(C)C(C)(C)O2)ccc1F

### **Products:**

1. CC(C)(C)c1cc(-c2ccnc(N)n2)ccc1F

Typical conditions: Pd catalyst.base.solvent

Protections: none

Yield: good

**Reference:** 10.1021/cr00039a007 and 10.1007/3418\_2012\_32 and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H

# 2.2.4 Reductive Amination of Aldehydes with Secondary Amines

### Substrates:

1. 6-Oxo-1,6-dihydropyridine-3-carbaldehyde - available at Sigma-Aldrich

2. 1-Ethylpiperazine - available at Sigma-Aldrich

# Products:

1. CCN1CCN(Cc2ccc(=O)[nH]c2)CC1

 $\textbf{Typical conditions:} \ \ \text{NaBH}(\text{OAc}) \\ 3 \ \ \text{or} \ \ \text{NaBH} \\ 3 \\ \text{CN}$ 

Protections: none

Yield: good

**Reference:** DOI: 10.1021/jo960057x and 10.1021/jm7009292 and 10.1073/pnas.1405685111 and 10.1002/ejoc.201101063 and 10.1038/ja.2017.61 and 10.1021/jm4013906

# 2.2.5 Bromination of aromatic compounds

### Substrates:

1. CC(C)(C)c1cc(-c2ccnc(N)n2)ccc1F

# Products:

1. CC(C)(C)c1cc(-c2nc(N)ncc2Br)ccc1F

Typical conditions: Br2.Fe

Protections: none

 $\bf Yield: \ good$ 

Reference: 10.1021/acs.accounts.6b00120

Retrosynthesis ID: 7777000

# 2.2.6 Amination of pyridones

### Substrates:

- 1. CC(C)(C)c1cc(-c2nc(N)ncc2Br)ccc1F
- $2. \ CCN1CCN(Cc2ccc(=O)[nH]c2)CC1 \\$

# Products:

 $1. \ CCN1CCN(Cc2ccc(Nc3ncc(Br)c(-c4ccc(F)c(C(C)(C)C)c4)n3)nc2)CC1 \\$ 

Typical conditions: 1.PCl5.2.amine

Protections: none
Yield: moderate

**Reference:** 10.1021/jm300780p AND 10.3390/molecules170910902 AND

10.1021/jm00392a017